



CARR Australia Pty Ltd

1 Spireton Place, Pendle Hill.
 PO Box 335, Pendle Hill.
 NSW 2145, Australia.
 Telephone 61-2-9636 9525
 Facsimile 61-2-9636 8844
www.carrgroup.com.au

CLEANING PRODUCTS APPROVED FOR USE ON CARRFLEX STANDARD

PRODUCT	CONCENTRATION
ACTIMAX	ACCORDING TO MANUFACTURERS SPEC
ALCO-DC-7	ACCORDING TO MANUFACTURERS SPEC
ALCO SANNI	ACCORDING TO MANUFACTURERS SPEC
ARJO DISINFECTANT	ACCORDING TO MANUFACTURERS SPEC
ASEPTI ACTIVE WIPES	READY TO USE
ASEPTIC ACTIVE LIQUID	ACCORDING TO MANUFACTURERS SPEC
AZO MAX MULTI SURFACE WIPES	ACCORDING TO MANUFACTURERS SPEC
AZO WIPE HARD SURFACE WIPES	ACCORDING TO MANUFACTURERS SPEC
CIDAL DISINFECTANT CLEANER	ACCORDING TO MANUFACTURERS SPEC
CLINELL HAND & SURFACE SANITISER DISINFECTANT	ACCORDING TO MANUFACTURERS SPEC
CLINELL UNIVERSAL SANITISING WIPES	ACCORDING TO MANUFACTURERS SPEC
CLINELL DETERGENT WIPES	ACCORDING TO MANUFACTURERS SPEC
CLINELL SPORICIDAL WIPES	ACCORDING TO MANUFACTURERS SPEC
DISINFECTING WIPES	ACCORDING TO MANUFACTURERS SPEC
DISPATCH DISINFECTANT	ACCORDING TO MANUFACTURERS SPEC
DIVERCLEANSE	125ML:4 LITRES WATER
ECOPINE	5% ECOPINE (1:20 WATER)
ISOWIPES (KIMBERLY-CLARK)	READY TO USE
JUMBO LIV-WIPES	ACCORDING TO MANUFACTURERS SPEC
LEGGESIDE HB (LEGGE SYSTEMS)	ACCORDING TO MANUFACTURERS SPEC
LEMEX	ACCORDING TO MANUFACTURERS SPEC
MEDIWIPES DISINFECTION (SULCO LIMITED)	ACCORDING TO MANUFACTURERS SPEC
MULTIKLEEN	ACCORDING TO MANUFACTURERS SPEC
NEUTRAL DETERGENT WIPES (REYNARD HEALTH SUPPLIES)	ACCORDING TO MANUFACTURERS SPEC
NON ACID BATHROOM CLEANER (CERTO)	ACCORDING TO MANUFACTURERS SPEC
OASIS PRO 12 CLEANER	5% (1 : 20 WATER)
ODACON SPRAY	ACCORDING TO MANUFACTURERS SPEC
PAL DETERGENT WIPES	ACCORDING TO MANUFACTURERS SPEC
PAL DISINFECTANT WIPES (GLOBALMEDICS)	ACCORDING TO MANUFACTURERS SPEC
PUREGREEN 24 (PURE GREEN)	ACCORDING TO MANUFACTURERS SPEC
ROOM CARE R2L	ACCORDING TO MANUFACTURERS SPEC
SPRING CLEAN (ECO LAB)	5% SPRING CLEAN (1:20 WATER)
STERRIMATT (PURE) FOAM	READY TO USE

SUPER SANI CLOTH- GERMICIDAL DISPOSABLE WIPES (PDI)	ACCORDING TO MANUFACTURERS SPEC
SURFACE SANITISING WIPES	ACCORDING TO MANUFACTURERS SPEC
TRIGENE (BIO STRATEGY)	ACCORDING TO MANUFACTURERS SPEC
TRIX (SOAP)	READY TO USE
TUFFIE 5 CLEANING WIPES	ACCORDING TO MANUFACTURERS SPEC
TURBO CLEANER (SUREHANDS)	ACCORDING TO MANUFACTURERS SPEC
USEALL (TRUE BLUE CHEMICALS)	ACCORDING TO MANUFACTURERS SPEC
VIRACLEAN (WHITELY MEDICAL)	ACCORDING TO MANUFACTURERS SPEC
VIREX (DIVERSEY)	ACCORDING TO MANUFACTURERS SPEC
VIRKON (DUPONT)	5% VIRKON (1:20 WATER)

IMPORTANT FOR ALL CLEANING PRODUCTS:

Wipe off cleaning product residue with a clean, soft cloth and rinse with fresh warm water immediately after use.

If particular cleaner/disinfectant is not listed we advise to test the particular surface disinfectant in advance or to contact us in case of any doubts.

Please contact Carr Australia Pty Ltd for further information.